

# Superior Testing at an Affordable Price

## 3300 Metals Testing Package

The Instron® Model 3380 100 kN Metals Testing System packages are an extremely affordable mechanical testing systems that contains all of the components necessary to start testing right away. Because of their ease of use and reliability, the 3300 Series is the ideal system for quality control and production environments. Systems are configured with Bluehill® Materials Testing Software and testing can easily start in two steps - open a method and press start.

### Key Features

- Method driven for simple or sophisticated tests
- World-class precision and accuracy
- Excellent long-term investment
- Safe and reliable
- Leasing options available
- Quick delivery

### Bluehill Software

Operators can setup a test using one of the pre-configured test method templates, gather results using the extensive calculations library, and generate an email report when the test is complete.

To begin a test simply open a method and press start – it's that easy. With multiple operators, varying skill levels, and continuous training needs, you want a user interface that is simple to learn. Bluehill provides all the capabilities you need to handle basic test requirements quickly and efficiently. Test control, data acquisition, plotting, calculations, and reporting are performed automatically. To further assist the operator, the Bluehill prompted method feature guides operators through repetitive test procedures with step-by-step instructions for all user levels.

For information on all packages available, visit [go.instron.com/3300packages](http://go.instron.com/3300packages) or call us at 1.800.564.8378



Metals Testing Packages include Bluehill Software, PC, Installation, and Training. Systems can be **configured** with a variety of grips, flex fixtures, and extensometers **to meet ASTM E8, ISO 6892, and more.**